



TIVA Administration Sets



BD Alaris[™] nexUS

PK

INFUSING

Ce	2.0ng/ml	Renifentanil Eleveld 10.0ng/ml
Cet	2.0ng/ml	31.7ml/h 7.0
		0.07ng/kg/min
		Cet 2.0ng/ml
		Cet 0.1h 23s

+ ADJUST -

MORE

BD Alaris[™] nexUS

PK

Propofol Eleveld 10.0ng/ml

Ce	1.8ng/ml	0.0ml/h 10.0
Cet	2.0ng/ml	0.00ng/kg/h
		Cet 4.0ng/ml
		Cet 0.3h 51s

+ ADJUST -

MORE

BD Alaris[™] nexUS

CC

INFUSING

15.0ml/h	BD PLASTIPAK 50	PRESSURE
0.25ng/min	VOLUME 0.56ml	1h 22m 38s
	0.56ng	TS=0.0
	VOLUME	VTBI

+ ADJUST -

FUROsemide (1 mg/ml)

BD Alaris[™] nexUS

Saline 0.9%

RATE 25.0 ml/h

VTBI 500 ml

0.6ml

19h 58m 21s

BD Alaris[™] nexUS

URGONocin

RATE 6.0ml/h

10.0ml/h

250 ml

0.2ml



24h

The full range of TIVA sets from BD.

500-002V

Qty per UOM: 50

- Set length 15cm
- 2-way Y connector
- Rotating male luer
- Polyethylene lined
- 1 pinch clamp



1

Gravity

Priming volume: 0.5ml

Tube inner diameter: 1.85mm

Tube outer diameter: 3.6mm

2

White

Priming volume: 0.6ml



Tube inner diameter: 1.7mm

Tube outer diameter: 3mm

500-003V

Qty per UOM: 50

- Set length 15cm
- 3-way Y connector
- Rotating male luer
- Polyethylene lined
- 2 pinch clamp



1

Gravity

Priming volume: 0.5ml

Tube inner diameter: 1.85mm

Tube outer diameter: 3.6mm

2

White

Priming volume: 0.6ml

Tube inner diameter: 1.7mm

Tube outer diameter: 3mm

3

Green

Priming volume: 0.3ml

Tube inner diameter: 0.5mm

Tube outer diameter: 2.5mm

Key to symbols



Female luer back-check valve



Female luer anti-siphon valve



Male luer

500-012V

Qty per UOM: 50

- Set length 220cm
- 2-way Y connector
- Rotating male luer
- Polyethylene lined
- 1 pinch clamp




1



1



1



1

Gravity

Priming volume: 0.5ml

Tube inner diameter: 1.85mm

Tube outer diameter: 3.6mm

2

White

Priming volume: 5.1ml

Tube inner diameter: 1.7mm

Tube outer diameter: 3mm

500-013V

Qty per UOM: 50

- Set length 190cm
- 3-way Y connector
- Rotating male luer
- Polyethylene lined
- 2 pinch clamp




1



2



1



1

Gravity

Priming volume: 0.5ml

Tube inner diameter: 1.85mm

Tube outer diameter: 3.6mm

2

White

Priming volume: 5ml

Tube inner diameter: 1.7mm

Tube outer diameter: 3mm

3


Green

Priming volume: 4.7ml

Tube inner diameter: 0.5mm


Tube outer diameter: 2.5mm

Product Code: 500-003V



- 1 Integrated anti-siphon valve¹
- 2 Colour-coded anti-siphon valves and clamps¹
- 3 Back-check valve¹
- 4 Multiple lumen hub¹
- 5 Polyethylene (PE) - lined tubing¹

Polyvinyl chloride (PVC) - free in the fluid path (no phthalates).



Reference:

1. BD technical data sheet: Extension sets for anaesthesia.

Additional safety information for TIVA sets, including Contraindications, Warnings, and Precautions may be found in the Instruction For Use (IFU).

Becton Dickinson Pty Ltd, Australia, Toll free: 1800 656 100
Becton Dickinson Limited, New Zealand, Toll free: 0800 572 468

bd.com

