



BD Totalys™ SlidePrep

Efficiency and confidence in slide processing

- Processes BD SurePath™ Liquid-based Pap test and Non-Gynaecological samples
- Integrates slide creation and staining from an enriched cell pellet
- Allows for random sample loading
- Fully integrates with the BD Totalys™ MultiProcessor to maintain continuous chain-of-custody



Help improve your lab's productivity

A new level of confidence through automation

- Consistent results with staining in individual chambers using fresh reagents for every sample
- Accurate reagent pipetting and precise timing for each sample reduces process variability
- Automatic barcode scanning ensures sample match and chain of custody

Efficient throughput

- Complete a single run of 48 slides in less than 62 minutes¹
- High throughput capability: up to 336 stained slides produced from cell pellet in an 8-hour shift²
- Minimal hands-on time allows staff to perform other tasks



BD Totalys™ SlidePrep

Set the precedent

Contact your BD representative today.

Specimen handling	
Specimen type	Enriched cell pellets and non-gynaecological cell pellets from the BD SurePath™ collection vial
On-board sample capacity	Up to 48 samples
Positive sample identification	2D barcodes on all centrifugation tubes and cytology slides ensure chain of custody throughout processing

Dimensions	
Height	86.4 cm (34 in)
Width	117 cm (46 in)
Depth	62.2 cm (24.5 in)
Weight	124.7 kg (275 lbs)
Separate bulk reagent station	35.6 cm H x 55.9 cm W x 35.6 cm D (14 in H x 22 in W x 14 in D)
Separate waste station dimensions	40.6 cm H x 33 cm W x 58.4 cm D (16 in H x 13 in W x 23 in D)

Environmental requirements	
Temperature	18–30 °C (64–86 °F)
Humidity	20–80% RH, non condensing
Altitude	Up to 2,000 m (6,562 ft)

Electrical requirements	
Input voltage	100–240 VAC
Rated current	12 A
Input line frequency	50–60 Hz
Power	1,000 Watts

1. BD Totalys™ SlidePrep Instrument User's Manual [500005570]. Franklin Lakes, NJ: Becton, Dickinson and Company; 2019.
2. Internal BD Document, PD-08129-67.

BD - Europe, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland

bd.com

